

Press Release


Super-photostable fluorescent labeling agent for super-resolution microscopy

~ A powerful tool for 3D and multicolor STED imaging of cellular ultrastructure ~

August 17, 2017

Chemists at ITbM, Nagoya University have developed a super-photostable fluorescent dye, PhoxBright 430 (PB430), to visualize cellular ultrastructure by super resolution microscopy. The exceptional photostability of this new dye enables continuous STED imaging and together with its ability to fluorescently label proteins, PB430 demonstrates its use in the 3D construction and multicolor imaging of biological structures.

PhoxBright 430 (PB430)




Molecular structure of the super-photostable fluorescent labeling agent, PhoxBright430 (PB430). PB430 consists of a fully ring-fused π -conjugated skeleton with an electron accepting phosphole *P*-oxide unit, which is responsible for its high photostability. The carboxylic acid of the molecule can be converted to the *N*-hydroxylsuccinimidyl (NHS) ester, which allows conjugation to biomolecules. Sulfonate groups ensure high solubility of the molecule in water.

Nagoya, Japan – Super resolution fluorescence microscopy, which has been acknowledged by the 2014 Nobel Prize in Chemistry, allows us to visualize biological systems and gain a detailed understanding of the complex dynamics of biomolecules. In particular, **stimulated emission depletion** (STED) microscopy is widely utilized to examine various processes in living systems, due to its fast acquisition speed and compatibility for various biological samples.


Chemists at the Institute of Transformative Bio-Molecules (ITbM) of Nagoya University have developed a new photostable fluorescent dye, PhoxBright 430 (PB430), which enables continuous super resolution STED imaging of fluorescently labeled cells. As PB430 contains a functional moiety to allow conjugation with an antibody, specific biomolecule targets within a cell can be stained by immunofluorescence. The exceptional photostability of PB430 has proven to be useful in constructing a 3D-STED image of microtubules (tube-shaped protein structures that form the cytoskeleton, which is the internal framework of a cell) and achieve multicolor STED imaging of fluorescently immunolabeled cytoskeletons by the combination of photostable PB430 and commercially available dyes.

The unique properties of PB430 make it a powerful tool to reveal the structures and functions of cells, and has the potential to be applied for prolonged visualization of the movement of organelles and molecules within cells. The results of this study are recently reported in the *Journal of the American Chemical Society*.


Super photostable fluorescent labeling reagent
PhoxBright 430 (PB430-NHS)




★ Exceptional photostability



★ Repeated imaging in the Z axis



★ Construction of super resolution 3D images





Possible to achieve super resolution (~150 nm) in the Z axis

Properties of the photostable dye, PB430. PB430 shows exceptional photostability compared to conventional dyes, as shown in the number of confocal images that can be taken without losing its relative fluorescence intensity. The new dye can also label proteins, such as antibodies and is able to scan the cell along the z-axis with a step of 50 nm, resulting in the re-construction of the 3D-STED image of the cytoskeleton with z-axis resolution of 160 nm.

In order to achieve high resolution in STED microscopy, a biological sample labeled with a fluorophore (a fluorescent dye used to label biological samples such as proteins, tissues and cells) is irradiated with a fluorescence excitation beam along with a doughnut-shaped STED beam to suppress the excitation of the surrounding molecules. Despite being a powerful technique, the need for a high intensity STED beam, which leads to photobleaching of dyes, has been a major issue for the practical use of STED microscopy in continuous live cell imaging.

C-Naphox



Molecular structure of the photostable fluorescent dye, C-Naphox. C-Naphox consists of a naphthophosphole P-oxide and a trimethylamine moiety. Carbon atoms are colored in gray, oxygen atoms in red, and a nitrogen atom in blue.

Researchers have previously reported a fluorescent dye, “C-Naphox”, which has exhibited strong photostability in STED imaging. The group has now optimized the structure of C-Naphox and developed a new photostable fluorescent dye PB430, which can be used for the visualization of structures within cells.

“C-Naphox has demonstrated to be an extremely photostable dye for the STED imaging of various materials, so we decided to tune its properties further so we can apply it to visualize ultrastructures,” says Masayasu Taki, an associate professor at ITbM and one of the leaders of this research. “Our new dye, PB430 shows increased solubility in water, fluoresces efficiently in aqueous media, and is able to label antibodies. We were pleased to observe that it also exhibits high photostability under STED conditions,” describes Taki.

This article “Super-Photostable Phosphole-Based Dye for Multiple-Acquisition Stimulated Emission Depletion Imaging” by Chenguang Wang, Masayasu Taki, Yoshikatsu Sato, Aiko Fukazawa, Tetsuya Higashiyama, and Shigehiro Yamaguchi is published online in the *Journal of the American Chemical Society*. DOI: 10.1021/jacs.7b04418 (<http://dx.doi.org/10.1021/jacs.7b04418>)

Author Contact

Professor Shigehiro Yamaguchi
Institute of Transformative Bio-Molecules (WPI-ITbM), Nagoya University
Furo-Cho, Chikusa-ku, Nagoya 464-8601, Japan
TEL: +81-52-789-2291
E-mail: yamaguchi@chem.nagoya-u.ac.jp
Homepage: <http://orgreact.chem.nagoya-u.ac.jp/Home.html>

Associate Professor Masayasu Taki
Institute of Transformative Bio-Molecules (WPI-ITbM), Nagoya University
Furo-Cho, Chikusa-ku, Nagoya 464-8601, Japan
TEL: +81-52-789-5248
E-mail: taki@itbm.nagoya-u.ac.jp

Media Contact

Dr. Ayako Miyazaki
Institute of Transformative Bio-Molecules (WPI-ITbM), Nagoya University
Furo-Cho, Chikusa-ku, Nagoya 464-8601, Japan
TEL: +81-52-789-4999 FAX: +81-52-789-3053
E-mail: press@itbm.nagoya-u.ac.jp

About WPI-ITbM (<http://www.itbm.nagoya-u.ac.jp/>)

The Institute of Transformative Bio-Molecules (ITbM) at Nagoya University in Japan is committed to advance the integration of synthetic chemistry, plant/animal biology and theoretical science, all of which are traditionally strong fields in the university. ITbM is one of the research centers of the Japanese MEXT (Ministry of Education, Culture, Sports, Science and Technology) program, the World Premier International Research Center Initiative (WPI). The aim of ITbM is to develop transformative bio-molecules, innovative functional molecules capable of bringing about fundamental change to biological science and technology. Research at ITbM is carried out in a "Mix Lab" style, where international young researchers from various fields work together side-by-side in the same lab, enabling interdisciplinary interaction. Through these endeavors, ITbM will create "transformative bio-molecules" that will dramatically change the way of research in chemistry, biology and other related fields to solve urgent problems, such as environmental issues, food production and medical technology that have a significant impact on the society.



From the left: Tetsuya Higashiyama, Aiko Fukazawa, Yoshikatsu Sato, Masayasu Taki, Chenguang Wang, Shigehiro Yamaguchi